

 <p><b>TARUMANAGARA</b> HUMAN CELL TECHNOLOGY LABORATORY THCT Lab</p>	<p><b>TARUMANAGARA HUMAN CELL TECHNOLOGY (THCT) LABORATORY</b>  <b>FACULTY OF MEDICINE, TARUMANAGARA UNIVERSITY</b>  <b>Jl. Letjen S. Parman No. 1; Jakarta 11440 INDONESIA</b>  Phone. +62 215696 3254, Fax. +62 21 5696 7325, Email. thctlab11@gmail.com</p>	
	<p><b>STANDARD OPERATING PROCEDURE</b> <b>(SOP)</b></p>	Document Num.: S-021
		Total pages: 3
	<p>Fibroblast Matrix Implant (FMI) for Hernia Repair</p>	Release date: 16 September 2019
Effective date: 16 September 2019		
Version: 1		
	First draft:	

### ***Scaffold fabrication***

- PLLA pellets (#B6002-2) are purchased from Durect Lactel (Cupertino, California, USA)
- Grind and Sieves sodium chloride (NaCl) particulates ranging in diameters from 355 to 425  $\mu\text{m}$
- Put 9 gr of the NaCl particulates into an aluminum pan (diameter of 40 mm)
- Dissolve 1 gr of PLLA pellets in 5.3 mL of chloroform & vortex the PLLA solutions
- Pour the PLLA solutions into the aluminum pans and mix with NaCl particulates
- Air-dry the mixtures in a fume hood at room temperature for more than 24 hours to evaporate the chloroform
- Detach the PLLA/NaCl composite from the aluminum pans
- Place the composite in a vacuum drying service (Yamato, Tokyo, Japan) at -0.1 mPa for another 3 days
- Wash the PLLA/NaCl composites in deionized water to leach out the NaCl particulates and the water is changed every hour for 20 times
- The PLLA scaffold (matrix) is formed after drying
- Cut the scaffold into discs with 1 cm diameter and 2 mm thickness for further experiments.

### ***Surface coating with collagen***

- Immerse the PLLA scaffold discs in the 1% (w/v) type I porcine collagen solution (Wako, Osaka, Japan) and deaerate under a reduced pressure (20 Pa) with a vacuum freeze dryer FDU-2200 (EYELA, Tokyo, Japan)
- Place the collagen/PLLA constructs onto 100  $\mu\text{m}$  cell-strainer membranes (Falcon, Durham, USA)
- Centrifuge at 2000 x g for 10 minutes at 4°C to remove collagen solution excess from the pores
- Freeze the collagen solution-coated scaffold discs at -80°C for a minimum of 4 hours
- Then, freeze-dry the scaffold discs under vacuum (<5 Pa for more than 24 hours)

### ***Plasma surface treatment and sterilization***

- Treat the PLLA collagen-coated discs with plasma for surface cleaning and activation with O<sub>2</sub> for 10 min (custom-made plasma treatment system, Diener electronic GmbH, Ebhausen, Germany)
- Seal the matrices in Tyvek® bag, and sterilize in the Diener machine with 30% H<sub>2</sub>O<sub>2</sub> solution (Merck, Darmstadt, Germany) at 30°C for 12 hours
- Then transfer the matrices into 24-well plates and store under vacuum at 4°C until use.

### ***Isolation of fibroblast cells from rats***

- The primary rat fibroblasts are isolated from abdominal skin specimens of approximately 4x3 cm Sprague Dawley rats (obtained from National Agency of Drug and Food Control, Jakarta)
- Under aseptic conditions, dissect and remove the abdominal skin
- Wash the tissue 5 times with PBS (Gibco, Grand Island, NY) containing 2% antibiotic-antimycotic (Sigma)
- Wash 5 times with EGTA solution (Sigma) containing 2% antibiotic-antimycotic
- Digest it for 5 hours with 30 mL of 200 U/mL collagenase type I (Gibco) solution at 37°C, with 5% CO<sub>2</sub>, and 95% relative humidity
- Gently scrap the dermis side of the skin with cell scraper
- Sieve the cells (1 mm and 100 µm mesh size)
- Wash twice and resuspend with Dulbecco's Modified Eagles Medium (Sigma, St.Louis, Missouri, USA) buffered with N-(2-hydroxyethyl)piperazine-N0-2-ethanesulfonic acid (HEPES) (Sigma), supplemented with 20% fetal bovine serum (Gibco, Grand Island, New York, USA) and 1% of antibiotic/antimycotic (Sigma)
- Culture the cells in a T75 cell culture flask
- Incubate the cells in incubator equilibrated with 5% CO<sub>2</sub> at 37°C
- When the cells have reached confluence, dissociate the cells with 2 ml 0.05% Trypsin-EDTA (Gibco) for 5 min in incubator equilibrated with 5% CO<sub>2</sub> at 37°C
- Add 5 ml of complete culture medium to stop the dissociation
- Centrifuge at 130 G for 5 min to pellet the cells
- Aspirate the medium and resuspend it with fresh culture medium
- Count the cells using the Tryphan blue exclusion method.

### ***Cytocompatibility tests***

- In order to test the scaffolds cytocompatibility, the PLLA collagen-coated matrices are seeded with primary rat dermal fibroblasts and cells viability is analyzed
- Seed the primary rat dermal fibroblasts onto the PLLA collagen coated matrices (500,000 fibroblasts/matrix)
- As an internal control, seed a same number of cells (500,000 fibroblasts) directly into tissue culture wells-plate
- Incubate at 37°C, with 5% CO<sub>2</sub> and 95% relative humidity for 24 hours
- Analyze the cell viability using Cell Counting Kit-8® viability assay (CCK-8) (Sigma)
- Measure the absorbance using the Multiskan reader at λ450 nm (Multiskan Ex, Thermo Scientific, USA)
- Results were expressed as fold change relative to control percent adhesion (cells adhered on the wells, 100%).

### ***Cells labelling and seeding***

- Label the fibroblasts with CellVue® NIR815 dye (Li-Cor #929-90020, Lincoln, Nebraska, USA) according to the kit manual before seeded onto the PLLA matrices
- Seed 100 µL of the fibroblast cell suspensions onto the PLLA scaffolds (500,000 fibroblasts/matrix)

- Add 300  $\mu$ L culture medium into each well, then after 1 hour add another 700  $\mu$ L culture medium into each well
- Incubate it in the incubator with 5% CO<sub>2</sub> at 37°C for overnight before implantation in order to allow the cells adhere to the scaffolds

***FMI in hernia animal model***

- Male Sprague-Dawley (SD) rats are used as hernia model by making 2 cm incision at the midline abdomen
- Separate peritoneum from the muscle to create space or pouch for the matrices implantation (both side of the midline)
- Close the peritoneum using 6-0 suture
- Implant or place the engineered fibroblast-matrices (2 pcs) at the abdominal area, between the abdominal muscle and peritoneum (sublay hernia repair method), one on each side.
- Close the muscle and the abdominal skin using 4-0 suture
- Close the surgical wound with tegaderm and hypafix
- Scan the implantation sites using Imaging machine (Li-Cor Pearl Trilogy Digital Imaging System) to observe the viability and distribution of fibroblast in the matrix at the implantation site
- After the end-point, anesthetize and euthanized the rats. Then, necropsy is performed to collect the tissue at implantation sites; store it inside the formalin solution for histopathological analysis

***Histopathology analysis***

- Fixate the specimens in 10% buffered formalin and prepare the histological slides
- Stain the histological slides with Hematoxylin and Eosin, Masson Trichome to observe the collagen formation, and Immunohistochemistry for fibroblast

Validation	Prepared by:	Checked by:	Authorized by:
Name	Jennifer Lheman	Dr. dr. Siufui Hendrawan, M.Biomed	Prof. Dr. Med. Hans U Baer
Signature			
Date	16 September 2019	16 September 2019	



REPUBLIK INDONESIA  
KEMENTERIAN HUKUM DAN HAK ASASI MANUSIA

# SURAT PENCATATAN CIPTAAN

Dalam rangka perlindungan ciptaan di bidang ilmu pengetahuan, seni dan sastra berdasarkan Undang-Undang Nomor 28 Tahun 2014 tentang Hak Cipta, dengan ini menerangkan:

Nomor dan tanggal permohonan : EC00201972907, 26 September 2019

**Pencipta**

Nama : **Dr. dr. Siufui Hendrawan, M.Biomed**

Alamat : Jl. Letjen S. Parman No.1, Universitas Tarumanagara, Fakultas Kedokteran, Jakarta Barat, DKI Jakarta, 11440

Kewarganegaraan : Indonesia

**Pemegang Hak Cipta**

Nama : **Dr. dr. Siufui Hendrawan, M.Biomed**

Alamat : Jl. Letjen S. Parman No.1, Universitas Tarumanagara, Fakultas Kedokteran, Jakarta Barat, DKI Jakarta, 11440

Kewarganegaraan : Indonesia

Jenis Ciptaan : **Buku Panduan/Petunjuk**

Judul Ciptaan : **Fibroblast Matrix Implant (FMI) For Hernia Repair**

Tanggal dan tempat diumumkan untuk pertama kali di wilayah Indonesia atau di luar wilayah Indonesia : 16 September 2019, di DKI Jakarta

Jangka waktu perlindungan : Berlaku selama hidup Pencipta dan terus berlangsung selama 70 (tujuh puluh) tahun setelah Pencipta meninggal dunia, dihitung mulai tanggal 1 Januari tahun berikutnya.

Nomor pencatatan : 000155958

adalah benar berdasarkan keterangan yang diberikan oleh Pemohon.  
Surat Pencatatan Hak Cipta atau produk Hak terkait ini sesuai dengan Pasal 72 Undang-Undang Nomor 28 Tahun 2014 tentang Hak Cipta.



a.n. MENTERI HUKUM DAN HAK ASASI MANUSIA  
DIREKTUR JENDERAL KEKAYAAN INTELEKTUAL

Dr. Freddy Harris, S.H., LL.M., ACCS.  
NIP. 196611181994031001